



Welch Allyn
Braun ThermoScan[®]
PRO 6000 Ear Thermometer

Advanced measurement technology that captures accurate, reliable temperature readings in a comfortable way for patients

WelchAllyn[®]

pharmaniaga[®]
passion for patients



Why measure in the ear?

Clinical studies have shown that the ear is an excellent site for measurement as temperatures taken in the ear reflect the body's core temperature ¹

Advantages of measuring the ear include ease of access and speed

- > Less invasive for the patient than oral, axillary or rectal temperature measurements
- > No mucus membrane contact
- > Staff can quickly capture measurements in the ear, even with challenging patient populations

¹Guyton A C, Textbook of medical physiology, W.B. Saunders, Philadelphia, 1996, p 919



The Welch Allyn Braun ThermoScan PRO 6000 Ear Thermometer has the speed and simplicity you need for accurate, reliable temperature readings.

Designed to improve staff productivity and patient satisfaction



Advanced Technology

- > Innovative PerfectTemp® technology adjusts for variability in probe placement
- > ExacTemp® technology detects stability of the probe during measurement
- > A pre-warmed sensor helps support accurate measurements
- > Designed with plastics compatible with common medical-grade cleaning products



Enhance Clinical Workflow

- > Memory recall button displays last measurement taken
- > 60-second pulse timer assists with manual measurement of pulse rate and respirations
- > C/F button for quick scale conversion of displayed reading
- > Electronic and mechanical security features help prevent theft and loss



Improve Patient Experience

- > Pre-warmed probe tip enhances patient comfort during reading
- > Probe covers are soft-walled for added comfort, and disposable to help minimize risk of cross-contamination
- > Temperature readings are quick and accurate—an ideal method for patients who are unable to cooperate during measurement



How does the Welch Allyn Braun ThermoScan PRO 6000 Ear Thermometer work?

When the Braun PRO 6000 is placed in the ear, it continuously monitors the infrared energy emitted by the tympanic membrane and surrounding tissues, until a temperature equilibrium has been reached and an accurate measurement can be taken.

PerfectTemp Technology— the Braun PRO 6000 advantage

PerfectTemp technology helps overcome challenges presented by ear canal anatomy and variability in probe placement. The thermometer collects information about the direction and depth of ear probe placement and automatically incorporates it into the temperature calculation. This helps support the accuracy of measurement as compared to core temperature, especially when probe positioning is not ideal.



PerfectTemp technology overcomes the potential for low readings, compared to core, by adjusting for factors that impact measurement accuracy.

Device options and accessories



Small Cradle (one-box)

Two cradles offer storage of probe cover boxes, and pop up for easy probe cover attachment. The large cradle can be wall mounted for added security when used with an optional security tether.



Large Cradle (two-box)



Optional Charging Station

Has storage for 200 probe covers, includes Rechargeable Battery Pack and enables electronic security settings when used with Welch Allyn Service Tool Software.



Optional Security Tether

Keeps thermometer attached to the cradle, to minimize loss and theft.

Select Specifications

Measurement Time	2–3 seconds
Proprietary Aiming Technology	PerfecTemp supports accuracy by detecting direction and depth of probe placement in ear and adjusts temperature calculation
	ExacTemp helps support reliability and accuracy by detecting the stability of probe placement during measurement
Displayed Temperature Range	68 °F to 108 °F (20.0 °C to 42.2 °C)
Operating Temperature Range	50 °F to 104 °F (10.0 °C to 40.0 °C)
Accuracy	±0.4 °F (±0.2 °C) for the range 95 °F to 107.6 °F (35.0 °C to 42.0 °C)
Display Type	Backlit liquid crystal display, four digits plus special icons
Automatic Power Off	Approximately 10 seconds after last measurement
Cradle Options	Small one-box (holds 20 covers), Large two-box (holds 40 covers)
Security Option—Electronic	Time interval countdown (requires Charging Station)
Security Option—Mechanical	Wall-mountable large cradle and optional security tethers
Other Features	Pulse timer, memory recall and C/F conversion buttons
Optional Accessories	Charging Station with Rechargeable Battery Pack, Security Tethers
Battery	2 x 1.5 V type AA (Alkaline) or Rechargeable Battery Pack
Battery Life	1000 measurements (AA Alkaline); 700 measurements on a fully charged Rechargeable Battery Pack (requires Charging Station); 350,000 measurements over the life of Rechargeable Battery Pack
Weight	3.6 oz (100 g) without batteries
Dimensions	6" × 1.7" × 1.3" (152 mm × 44 mm × 33 mm)
Standard	ASTM 1965-98, EN 12470-5:2003, MDD (93/42/EEC) Annex II
Warranty	3 years (Thermometer, Charging Station and Rechargeable Battery Pack)

Ordering Information

- 06000-200 Braun ThermoScan PRO 6000 and Small/One-box Cradle
 06000-300 Braun ThermoScan PRO 6000 and Large/Two-box Cradle

Accessories

- 06000-005 Braun ThermoScan PRO 6000 Disposable Probe Covers (5000/case)
 06000-100 Braun ThermoScan PRO 6000 Charging Station (includes Rechargeable Battery Pack)
 106201 Security Tether, 6 ft
 106204 Security Tether, 9 ft
 104894 Braun ThermoScan PRO 6000 Rechargeable Battery Pack
 01802-110 Model 9600 Plus Temperature Calibration Tester



WelchAllyn®

Authorised Malaysia Distributor
 Pharmaniaga Logistics Sdn Bhd
 No.7, Lorong Keluli 1B, Kawasan Perindustrian Bukit Raja Selatan,
 Seksyen 7, 40000 Shah Alam, Selangor Darul Ehsan,

General line : 03-3342 9999
 Customer Care Line : 1-800-888-313
 Fax : 03 -3342 5715



Welch Allyn

CareTemp® Touch Free Thermometer

A quick and convenient way to screen patient temperatures without contact



In fast-paced environments where speed and efficiency matter, the CareTemp Touch Free Thermometer is a valuable screening tool for many clinical settings.

Support a safe and comfortable patient experience

- Infrared sensor reads temperature 4-6 cm from center of forehead—no patient skin contact required
- Quickly obtain readings with minimal patient involvement to determine if a patient presents with a normal or elevated temperature
- Uses a touch free sensor head to help minimize the risk of cross-contamination

Enhance staff productivity and satisfaction

- Measurement feedback system detects proper operating distance to help provide consistent accuracy
- Device head does not require cleaning between each patient to help save time and cost
- Features to help enhance workflow include: pulse timer, temperature memory recall button and a backlit display

CareTemp® Touch Free Thermometer

Large, backlit LCD display is intuitive and easy to read

Dedicated memory recall button to display last measurement taken

Ergonomic handle design for a comfortable and secure hold



Proprietary distance sensor assists with consistent measurements

Green LED illuminates when correct operating distance has been achieved

Easily switch displayed reading between Celsius and Fahrenheit scales

60-second timer for manual pulse and respiration measurement



Ideal Care Environments*

- > Triage
- > Emergency Department
- > Pediatric Clinic
- > Urgent Care
- > School Health Services
- > Long-term Care
- > Mental Health

**For care environments when the forehead site may not track with core temperature changes, Welch Allyn recommends SureTemp® Plus (oral/axillary/rectal) and Braun ThermoScan® (ear) thermometers.*

Ordering Information

105801 CareTemp Touch Free Thermometer

Standards

ASTM 1965-98, ISO 80601-2-56, MDD (93/42/EEC) Annex II

Specifications

Device Type	Forehead touch free infrared thermometer
Measurement Distance	4 to 6 cm from the center of the forehead
Aiming Technology	Measurement feedback system detects when the proper operating distance is achieved to help provide consistent accuracy
Displayed Temperature Range	34 to 43 °C (93.2 to 109.4 °F)
Display Type	Backlit liquid crystal display, four digits plus special icons
Display Resolution	0.1 °C (0.1 °F)
Laboratory Accuracy	±0.2 °C (±0.4 °F) for the range 34 to 43 °C (93.2 to 109.4 °F)
Operating Temperature	16 to 40 °C (60.8 to 104 °F)
Operating Humidity	Up to 85% non condensing
Storage Temperature	-25 to 55 °C (-13 to 131 °F)
Storage Humidity	Up to 85% non condensing
Automatic Power Off	Approximately 10 seconds after last measurement
Battery	2 × 1.5 V AAA (LR03)
Battery Life	Approximately 3,300 measurements
Weight	78 g with batteries
Dimensions	159.3 mm x 41.3 mm x 55.8 mm
Warranty	2 years

Authorized Malaysia Distributor
Pharmaniaga Logistics Sdn Bhd
No.7, Lorong Keluli 1B, Kawasan Perindustrian Bukit Raja Selatan,
Seksyen 7, 40000 Shah Alam, Selangor Darul Ehsan,

General line : 03-3342 9999
Customer Care Line : 1-800-888-313
Fax : 03 -3342 5715



WelchAllyn®