

Sinovac Covid-19 1-Dose Vaccine Extraction Method



Syringe Size Requirement
Low Dead Volume (LDV)
1ml

Recommended Needle Size
**21G with the length at least
38mm (1 ½")**

Note: The needle length is important to ensure the needle tip could reach at bottom of the vial

Vial Position
Upright position at 45° tilted



Precaution During Withdrawal

1. Make sure the **correct syringe size** is used as usage of an oversized syringe will result in inaccurate content reading for small fill volume.
2. **Avoid over priming.** Over priming could cause the syringe's vaccine content to spill out, rendering the remaining vaccine content in the syringe to be below 0.5ml.
3. If there is insufficient volume withdrawn, check if there is any more vaccine residual leftover in the vial at upright position for at least 30 seconds.