


Sinovac Covid-19 2-Dose Vaccine Extraction Method



Syringe Size Requirement
Low Dead Volume (LDV)

1ml

Recommended Needle Size

Recommended Needle Size

21G with the length at least 38mm (1½")

Note: The needle length is important to ensure the needle tip could reach at bottom of the vial

Vial Position
Upright position at 45° tilted



Precaution During Withdrawal

1. For the first dose of vaccine, extract approximately 0.58-0.59ml of vaccine of vaccine from the vial.
2. The needle for extraction (21G) can be left in the vial to be used for the extraction of the 2nd dose of vaccine.
3. After the needle of the syringe is changed to 23G/25G with the length of 25mm and "prime" to the 0.50ml level on the syringe to start the injection.
4. For the 2nd dose of vaccine, extract all the remaining vaccine in the vial and repeat step-3 above.
5. If there is insufficient volume withdrawn for 2nd dose, check if there is any more vaccine residual leftover in the vial at upright position for at least 30 seconds.